

Mike's Moss Transformation Protocol

(Compiled from various sources, mostly Cove protocol)

Updated 2/11/2018

I. Preparation

- A. 8.5% D-Mannitol Solution (8+ mL/transformation + ~30 mL/protoplast strain)
- B. 10.28% D-Mannitol Solution (10 mL/protoplast strain)
- C. CaPW (8.5% D-Mannitol + 10 mM CaCl₂; autoclave) (10 mL/protoplast strain)
- D. 1% MES pH 5.6 (autoclave) (1 mL/protoplast strain)
- E. 1M MgCl₂ (autoclave) (<1 mL/protoplast strain)
- F. 1M Ca(NO₃)₂ (autoclave) (<1 mL/protoplast strain)
- G. 1M Tris•HCl pH 8 (autoclave) (50 µL/~20 transformations)
- H. Autoclave 2.0 g PEG-6000 in a 15 mL tube (1 tube/~20 transformations)
- I. PRMT medium BCD + 2.5 mM DAT, 8% D-mannitol, 0.4% Agar + 10 mM CaCl₂ (@50°C)
- J. PRMB plates BCD + 2.5 mM DAT, 6% D-mannitol, 0.8% Agar, 10 mM CaCl₂ (@50°C)
- K. 5-day-old moss chloronema grown on BCD+DAT + 0.6 mM CaCl₂ (± glucose)
- L. 10 to 30 µg linearized plasmid DNA in <30 µL (Digest, extract, ethanol precipitate, dry)
- M. BCD+DAT Medium (700 ml): 700 µl 500 mM MgSO₄; 700 µl 616 mM KH₂PO₄/303 mM K₂HPO₄ pH 6.5; 7 ml 1 M KNO₃; 0.7 ml Trace element solution; 8.75 mg FeSO₄•6H₂O; 5.6 g Agar; 0.644 g diammonium tartrate; 1 mM CaCl₂ (add @50°C)

II. Day 1 Before starting, set water baths to 45°C and 65°C, start DNA precipitation.

- A. Prepare 2% Driselase Solution in 8.5% D-Mannitol (0.2g/10ml) 15 minutes at RT, Centrifuge 2500 g, Filter sterilize supernatant.
 - B. Collect protonema to a sterile 50 ml tube, pipet away excess liquid, and add ~15 mL 8.5% D-Mannitol per plate of tissue.
 - C₁. Add 1/3 volumes of 2% Driselase Solution (0.5% final). Incubate 30 to 60 at RT with occasional gentle shaking (<2 hours). While waiting, complete steps C₂ to C₅.
 - C₂. Start PEGT preparation: add 4.45 mL 8.5% D-mannitol 500 µL Ca(NO₃)₂ 50 µL Tris pH 8, to a sterile PEG tube. Dissolve at 65° then allow to cool completely before use.
 - C₃. Pellet the plasmid DNA, rinse with 70% ethanol, aspirate, dry, and resuspend to 1 µg/µl
 - C₄. Melt PRMT and add 10mM CaCl₂
 - C₅. Label and place cellophanes on plates
 - D. Filter digested tissue through 70 µ mesh basket. (or 1st through 100µ basket then 40/70µ)
 - E. Pellet at 100-200 g for 4 minutes with no braking (1000 RPM in CL2 centrifuge).
 - F. Resuspend in 10 mL CaPW. Repeat once but use 8.5% D-Mannitol.
 - G. Count protoplasts and prepare enough fresh MMM Medium
Total # protoplasts = [# protoplasts in 1 mm² (5x5) grid] x [# mL protoplasts] x 10⁴
MMM = 8.85 mL 10.28% D-Mannitol, 150 µL 1M MgCl₂, 1 mL 1% MES pH 5.6.
 - H. Pellet as before and resuspend in MMM at 1.67x10⁶ cells/mL
 - I. Add <30 µL linearized plasmid DNA to 15 mL conical tubes (or round-bottom tubes)
 - J. Add 300 µL protoplasts
 - K. Add 300 µL PEGT (5 minutes)
 - L. 45°C for 5 minutes
 - M. 25°C for 10 minutes
 - N. Add 300 µL 8.5% D-Mannitol and swirl, repeat 4 additional times at ≥1 intervals.
 - O. Add 1 mL 8.5% D-Mannitol and swirl, repeat 4 additional times at ≥3 intervals.
 - P. Pellet as before (if round-bottom tubes were used, transfer to conical tubes.)
 - Q. Resuspend in 1 mL 8.5% D-mannitol.
 - R. Add 5 mL molten PRMT (45°C).
 - S. Dispense 2 mL per PRMB+cellophane plates.
 - T. Put in growth room
- III. Day 5 or 6: Transfer cellophane to Selective BCD(±DAT) plates (grow 10-17 days)
20-50 µg/ml G418, 15-30 µg/ml Hygromycin, 50 µg/ml Zeocin (replace weekly?), 100-150 µg/ml Gentamicin
- IV. Day ~20: Move (or pick resistant protonemata) to non-selective media (grow 7-14 days)
- V. Day ~34: Pick a few cells at 2-3 edges to selective media (grow 10+ days)
- VI. Day ~48: Check insertion site PCR (5' native/foreign, 3' native/foreign, empty site)