

Ultra Competent *E. coli* Cells

- Inoculate 250 mL **SOB** media with 10 - 12 colonies from a freshly streaked plate.
- Grow to an OD₆₀₀ of 0.94 at 18°C (or RT). (Doubling time is 4-10 hours). (1/04: 9.7 hr for INVaF')
- Chill on ice for about 10 minutes

****From this point on, it is critical to keep cells below 4°C****

****Work in cold room when possible****

- Pellet cells at 2500 x g for 10 minutes at 4°C.
- Resuspend cells gently in 80 mL of *ice cold TB*.
- Leave on ice for ~10 min.
- Pellet cells as before
- Resuspend cells gently in 10 mL *ice cold TB*.
- Add 10 mL *ice cold TB + 14% DMSO*.
- Wait 10 minutes, then aliquot into *pre-cooled* 1.5 mL tubes, and freeze them in liquid nitrogen.
- Store in -70 freezer.

Recipes

SOB+glycine

1 L

20 g Tryptone
5 g Yeast Extract
0.5 g NaCl
625 µL 4 M KCl
10 g glycine
pH to 7 with 1 M HCl (~800 µL)
Autoclave
10 mL 1 M MgCl₂

TB

100 mL

*Use the purest water available, and pure-water-rinsed glassware in preparation
0.30 g (10 mM) PIPES
0.22 g (15 mM) CaCl₂•2H₂O
1.86 g (250 mM) KCl
pH to 6.7 with 150-200(?) µL 10 M KOH
1.09 g MnCl₂•4H₂O
Filter sterilize

Things to autoclave

Centrifuge bottles
100 mL Graduated Cylinder (Chill)
1.5 or 2 mL eppendorf tubes (Chill)

Efficiency test: Transform 0.1 pg pUC19: 100 colonies ⇔ 10⁹ cfu/µg
1 pg pUC19: 100 colonies ⇔ 10⁸ cfu/µg
10 pg pUC19: 100 colonies ⇔ 10⁷ cfu/µg

Based upon Inoue, Nojima, and Okayama. 1990. *Gene* **96**: 23-28 with modifications suggested by Tang *et al.* 1994. *NAR* **22**: 2857-2858 and Akhtar, et al. 2000 *Anal. Biochem.* **277**: 273-276.

Transformation

- Thaw cells on ice. (If desired, split comp cells to pre-cooled tubes. I have gotten by with 12.5 μL /transformation.)
 - (Optional for ~2-3X efficiency increase) Add 4 μL of 0.5 **M** (or 0.14 μL of 100%) β -mercaptoethanol to cells per 100 μL competent cells. (20 **mM** final)
- Add DNA (I add <3 μL per 100 μL comp cells)
- Incubate in ice for ~30 minutes. (Longer incubations seem to help for plasmids >10 kb.)
- Heat shock at 42°C for 30 seconds then place on ice for another 2 minutes.
- Add 0.9 mL pre-warmed SOC media. (LB works but might affect efficiency)
- 37°C for 45 minutes (or longer). (Some people suggest shaking, others say it kills the cells)
- Plate.

From Inoue, Nojima, and Okayama. 1990. *Gene* **96**: 23-28;
 β -ME from Invitrogen's and Stratogene's protocols

“five minute” protocol

notes: This protocol is reportedly sensitive to salt in the DNA mix (50% less colonies at 5**mM** relative to 0**mM**). Although others claim to get twice the efficiency with this fast protocol, in my hands it was the same as that from above. For transforming supercoiled DNA, I squirt the transformation towards edge of plate and streak for single colonies. Do not use for kanamycin selection—only amp.

- Prewarm plates to 37°C for \geq 30 minutes.
- Thaw cells and add DNA.
- After >5 minutes, spread onto plate

From Pope & Kent *NAR* 24:536-7