

Mike's Moss Transformation Protocol

(Compiled from various sources, mostly Cove Lab protocol)


Updated Aug 7, 2024

I. Preparation (* = prep each time)

- A.* 10 to 30 μg plasmid DNA in $\leq 30 \mu\text{l}$.
I use Zymo midipreps for plasmids; it yields $200 \mu\text{l}$ $1.0 \pm 0.4 \mu\text{g}/\mu\text{l}$ plasmid. For oligo-mediated edits, use $10 \mu\text{l}$ of $25 \mu\text{M}$ double-stranded oligos after annealing with PCR machine. For transient/CRISPR transformations, use supercoiled plasmid.
- B.* 7-day-old moss chloronema grown on BCDAT (1 plate/7 transformations, at best)
The moss should be 4-7 days old (from last chop) and not too dense. The more chloronemata, the better. If the starting material had a lot of gametophores, you might need to go through two rounds of chopping.
- C.* Autoclave 2.0 g PEG-6000 in a 15 mL tube (1 tube/~20 transformations)
Add to 15 ml tube. Autoclave with cap loose covered with foil.
- D.* 8.5% D-Mannitol Solution (~15 mL/transformation + 50 ml/protoplast strain)
- E. 10.3% D-Mannitol Solution (9 ml/protoplast strain)
- F. 1% MES pH 5.6 (autoclave) (frozen aliquots in freezer) (1 ml/protoplast strain)
- G. 1 M MgCl_2 (autoclave) (frozen aliquots in freezer) (150 μl /protoplast strain)
- H. 1 M Tris•HCl pH 8 (autoclave) (frozen aliquots in freezer) (50 μl /~20 transformations)
- I. 1 M $\text{Ca}(\text{NO}_3)_2$ (autoclave) (500 μl /~20 transformations)
- J.* PRML medium=BCD + 2.5 mM DAT, 8% D-mannitol, 0.08% Agar, 10 mM CaCl_2 (@50°C)
or PRMT medium=BCD + 2.5 mM DAT, 8% D-mannitol, 0.5% Agar, 10 mM CaCl_2 (@50°C)
or 1.5% alginate+ 8.5% 8.5% D-mannitol plus 8.5% D-mannitol + 10 mM CaCl_2
- K.* PRMB plates BCD + 2.5 mM DAT, 6% D-mannitol, 0.8% Agar, 10 mM CaCl_2 (@50°C)
(need 3 or 4 plates per transformation plus 1-2 plates per starting strain.)
The more transformants you expect, the more plates you need. I generally split single-plasmid transformations to four plates and split two/three-plasmid transformations to three plates.
- L.* BCDAT selection plates (not needed until Day 5-8)

II. Day 1 Before starting, set water baths to 45°C and 65°C and turn on burner.

- A. Prepare 20 mg/ml Driselase + 3 mg/ml Cellulase R-10 Solution in 8.5% D-Mannitol.
Leave 15 minutes at RT, Centrifuge 2500 g for 10 min. Filter sterilize (0.22 μ) supernatant.
Take driselase out of freezer and warm it up before opening so that water doesn't condense inside the bottle. I make 5 ml if I only have a few transformations and 10 ml if I need more. Adding Cellulase R-10 seems to improve digestion. I leave them to dissolve in the hood with the light off, inverting the tube occasionally. I usually filter-sterilize the driselase/cellulase directly to the moss+mannitol tube in step C to save a tube.
- B. Collect protonema to a sterile 50 ml tube and add ~5 ml 8.5% D-Mannitol per plate of tissue using cell scraper and/or forceps.
Choosing how many plates of moss to harvest is tricky. If all goes well, you should get enough for at least 7 transformations per plate. I tend to start with more than enough if available, but this can be wasteful. I add 3 ml 8.5% Mannitol per ml of driselase solution being prepared.
- C₁. Add 1/3 volume of 20 mg/ml Driselase Solution (5 mg/ml final). Incubate 30 to 60 at RT with occasional gentle shaking (<2 hours). While waiting, complete steps C₂.
I set the tube in a clean tray on the shaker under the centrifuge set at the slowest speed. Protoplasts are fragile. If you don't see clumps breaking up within 30 minutes, there may be a problem.
- C₂. Start PEG/T preparation: add 4.45 ml 8.5% D-mannitol 500 μl $\text{Ca}(\text{NO}_3)_2$ 50 μl Tris pH 8, to a sterile PEG tube. Dissolve at 65° then allow to cool completely before use.
I combine the ingredients and put in the 65°C bath, eat lunch, then return and vortex the PEG/T. It should be mostly-but-not-completely melted/dissolved. Return to 65°C for ~10 more minutes and vortex until completely dissolved. Chill—it needs to be at room temperature by step L.
- D. Filter digested tissue through 70 μ strainer. (or 1st through 100 μ strainer then 40 μ)
I usually start the filtering the protoplasts after 1 h. I generally filter through a 70 μ strainer, but it sometimes clogs, taking a lot more time to filter. Filtering first with a 100 μ strainer, then through a 40 μ strainer seems to clog less and results in less debris in the prep (but wastes more plastic). Use heat-sterilized forceps to place the strainer(s) into new 50 ml Corning tubes, and add digested cells using a sterile transfer pipet.
- E. Pellet at 100-200 g for 4 minutes with no braking (1000 RPM in CL2 centrifuge = 167g).
I use the CL2 centrifuge at 1000 RPM. CJ used the Eppendorf at 250g.

- F. Pour off supernatant in one motion. Resuspend in 10 ml **8.5% D-Mannitol**.
Immediately after it stops, take tube out and pour into a waste beaker in one motion then recap immediately to limit contamination. If you hesitate, the pellet will be partially resuspended.
- G. Repeat centrifugation and resuspension. CJ did a 3rd wash.
- H. Count **protoplasts** and prepare enough fresh MMM Medium
Volume (ml) MMM needed = average hemocytometer count \div 17
The hemocytometer and counter are in a box above the Eppendorf centrifuge. Swirl tube to get protoplasts fully suspended, sterily add $\sim 20 \mu\text{l}$ to both halves of the hemocytometer, then cover with coverslip. Under the microscope, locate one of the central $1 \text{ mm} \times 1 \text{ mm}$ grids and count the intact (spherical) protoplasts. Repeat with the other $1 \text{ mm} \times 1 \text{ mm}$ grid. *Math: the coverslip sits $100 \mu\text{m}$ over the grid, so the volume ($1 \text{ mm} \times 1 \text{ mm} \times 0.1 \text{ mm}$) is $0.1 \mu\text{l}$. Multiply the average of the two counts by $10^4 \times 10$ to get the total number of protoplasts in the 10 ml. Divide that number by 1.7×10^6 to calculate the number of ml of MMM needed to get 1.7×10^6 protoplasts/ml.*
- 
- I. Pellet as before and resuspend in calculated volume of MMM.
10 ml MMM = 8.85 ml **10.3% D-Mannitol**, $150 \mu\text{l}$ 1 M MgCl_2 , 1 ml 1% MES pH 5.6.
- J. Add $\leq 30 \mu\text{l}$ plasmid DNA to 15 mL conical tubes, quick spin to get to tip
I pipet the DNA(s), protoplasts, and PEG to successively higher gradation marks on the tubes so that I don't have to change tips between samples. I use polystyrene tubes (Falcon 352095) instead of standard polypropylene tubes because cells may stick to them less.
- K. Add $300 \mu\text{l}$ **protoplasts**
- L. Add $300 \mu\text{l}$ **PEG/T** and mix gently
After adding the PEG/T to each sample, I gently swirl and tip the tube nearly horizontal to mix in droplets stuck to side of tube
- M. **45°C for 5 minutes**
Use smaller metal tube rack. While waiting, add $\sim 2''$ water to an ice bucket and add ice until 20°C
- N. **20°C for 5-10 minutes**
After ~ 5 minutes, I dry the tubes with paper towel and put back in a dry rack in the hood.
- O. Slowly add 10 ml **8.5% D-Mannitol** and mix by swirling.
I gently invert them until the cells in the tip become unstuck. Leave tubes horizontal rather than vertical.
- P. Incubate 1 hour (30 minutes to a few hours is okay). Meanwhile, label PRMB plates then add cellophanes.
Label them with at least the transformation number. I usually add a short description and the date.
- Q. While waiting, melt the PRMT if you are using it instead of PRML, then add CaCl_2 .
- R. Pellet as before. I do 4 at a time.
- S₁. If using PRMT (top agar), gently resuspend the pellets with $500 \mu\text{l}$ **8.5% D-Mannitol** and get the plates ready for the first transformation. Add [# plates – 0.5] ml molten PRMT to the tube and immediately suck it back into the pipet and drip 1ml across each plate.
I usually drip ~ 8 drops to different areas of each plate and then spread around with tip before it solidifies.
- S₂. If using alginate, gently resuspend the pellets with 1.5 or 2 ml of 10 mM CaCl_2 + 8.5% D-Mannitol (depending on whether you're plating to 3 or 4 plates). Add an equal volume of 1.5% alginate + 8.5% D-Mannitol to the cells then pipet 1 ml per plate.
The calcium causes the alginate to gel—add more CaCl_2 if you want it to gel faster/firmer, add less if it gels too quickly. With Sigma's A1112, 5 mM CaCl_2 gels well enough, but 3.75 mM did not. (Pierre-François Perroud uses for the alginate they use.) After moving to BCDAT plates (1 mM CaCl_2), the alginate seems to liquify some making cellophane transfers tricky. Don't use Fisher's 17777-2500; it gelled way too quickly with even 2.5 mM CaCl_2 and the gel disk shrunk $\sim 60\%$.
- S₃. If using PRML, resuspend with 2 ml (4 plates) or 1.5 ml (3 plates) and dispense $500 \mu\text{l}$ per plate.
PRML is faster and easier than PRMT but lowers the regeneration rate. I use it when I'm not worried about getting enough transformants. I use a 5 ml pipet to add 1.5 or 2 ml per transformation then use the same pipet to resuspend and pipet $500 \mu\text{l}$ to each plate. I use a transfer pipet or P1000 to pipet the others.
- T. *Optional* I also recommend plating $\sim 100 \mu\text{l}$ of the untransformed protoplast strain to a PRMB plate using the same method. This control is useful if you see poor regeneration or if the selection doesn't seem to have worked.
- U. Tape and put in growth chamber

III. Day 5-8: Transfer cellophane to BCDAT selection plates: *25 µg/ml G418, 15 µg/ml Hygromycin, 50 µg/ml Zeocin, 75-100 µg/ml Gentamicin, 125-150 µg/ml kanamycin.*

I always had moved cellophanes to selection after 4 days on PRMB, but Pierre-François' lab says you get more transformants if you wait 7 days. This increased transformation rate may be true, but I also saw an increased rate of plasmids meant to be maintained transiently that were stably integrated somewhere in the genome. For the no-DNA control plate, I cut the cellophane with a heat-sterilized scalpel and move strips to a plain BCDAT plate and to plates with each antibiotic. This has proven to be a helpful control a few times.

CRISPR deletions and CRISPR+oligo-mediated editing:

As soon as resistant plants are distinguishable (~week on selection), pick 24-96 to non-selective BCDAT. After growing for 7-14 days, pick a few cells from edge of each protonema to a new plate (reduces mosaicism). After growing for 7-14 more days, transfer a few cells to a new plate and use the rest for DNA preps for genotyping. PCR across the deletion/edit site. Confirm absence of the Cas9 plasmid by PCR or antibiotic resistance.

Targeted gene replacement with integrated selection cassette (with or without CRISPR facilitation):

After 7 days on selective media, transfer either the whole cellophanes or many transformants to non-selective media. After growing for 7-14 days, pick a few cells from edge of each protonema to a new selection plate (for the inserted resistance cassette, not the Cas9 plasmid's) and grow for 10+ days. Transfer a few cells to a new plate and use the rest for DNA preps for genotyping. PCR across both sides of the insertion site and confirm the absence of the deleted sequence and of Cas9.